

# DURABILITY II

## QUICK REFERENCE GUIDE



DURABILITY II is the first controlled study to focus on treating long, complex lesions with a single stent and to specifically test the performance of a 200 mm stent in the superficial femoral artery. After 36 months, DURABILITY II evidence supports a single stent strategy — implanting one long stent reduces fracture rates and increases vessel patency.

### STUDY OBJECTIVE

The study objective was to compare PTA and primary stenting using the EverFlex™ self-expanding peripheral stent system for the treatment of lesions in the native SFA or SFA and proximal popliteal arteries.

### STUDY DESIGN

- U.S. FDA IDE study
- Prospective, multicenter, non-randomized, single-arm study
- 287 subjects
- Independently adjudicated

### PRIMARY ENDPOINTS

- Safety:** Major adverse events at 30 days
- Effectiveness:** Primary patency at 1 year

### STUDY PARAMETERS

PARAMETER	Lesion Characteristics (n = 287)
Mean Lesion Length (core lab)	89.1 mm
Total Occlusions	48.1%
Calcification (moderate to severe)	70.0%
Moderate	26.8%
Severe	43.2%

### STUDY RESULTS

The 36-month results offer evidence that even in long, complex lesions, the EverFlex stent is able to sustain patency and durability.

DURABILITY II	RESULTS		
	12-mo <sup>1</sup>	24-mo <sup>2</sup>	36-mo <sup>3</sup>
Freedom from loss of primary patency (PSVR < 2.0) <sup>1</sup>	77.2%	66.0%	60.0%
Patency in lesions ≤ 80 mm	86.2%	80.8%	71.0%
Patency in lesions > 80 mm	69.6%	53.1%	50.5%
Fracture rate	0.4%	0.9%	0.9%

FREEDOM FROM TLR	12-mo <sup>1</sup>	24-mo <sup>2</sup>	36-mo <sup>3</sup>
All Subjects (n = 287)	86.0%	75.0%	70.0%
Lesions ≤ 80 mm (n = 133)	92.0%	89.0%	80.0%
Lesions > 80 mm (n = 154)	81.0%	64.0%	61.0%

<sup>1</sup>Kaplan-Meier analysis.

Confidential. For internal use only.  
Do not distribute.

**Medtronic**

# PERIPHERAL BARE METAL STENTS FOR FEMOROPOPLITEAL DISEASE

## CLINICAL DATA QUICK REFERENCE GUIDE<sup>†</sup>

	DURABILITY II <sup>1,2</sup>	RESILIENT <sup>3,4</sup>	STROLL <sup>5,6</sup>	SuperNOVA <sup>7</sup>	SUPERB <sup>8</sup>	OSPREY <sup>9,10</sup>
<b>Stent</b>	EverFlex™ Peripheral Stent	LifeStent™™ Vascular Stent	S.M.A.R.T.™™ Vascular Stent	Innova™™ Self-expanding Stent	Supera™™ Peripheral Stent	R2P™™ MISAGO™™ Stent
<b>Mean Lesion Length (mm)</b>	89.1	70.5	77.3	93.2	77.7	83.8
<b>Calcification (moderate to severe)</b>	70.0% (43.2% severe)	35.3% (24.5% severe)	19.3% (severe only)	Not Reported (35.6% severe)	72.8% (44.9% severe)	66.6% (31.4% severe)
<b>CTOs</b>	48.1%	17.0%	23.6%	44.7%	24.7%	Not Reported
<b>Patency</b>	77.2% (12 months) 66.0% (24 months) 60.0% (36 months)  Lesions ≤ 80 mm: 86.2% (12 months) 80.8% (24 months) 71.0% (36 months)	81.3% (12 months)  Lesions ≤ 80 mm: 85.0% (12 months)	87.9% (12 months) 80.5% (24 months) 72.7% (36 months)	66.4% (12 months)	86.1% (12 months)	82.9% (12 months)
<b>Fracture Rate</b>	0.4% (12 months) 0.9% (24 months) 0.9% (36 months)	3.1% (12 months) 4.1% (18 months)	2.0% (12 months) 2.3% (24 months) 3.6% (36 months)	1.9% (12 months) 2.2% (24 months)	0% (12 months) 0% (24 months) 0.6% (36 months)	0.5% (12 months) 2.0% (24 months) 2.0% (36 months)
<b>Key Message**</b>	Long, complex lesions, PSVR < 2.0	Shorter, less complex lesions, PSVR < 2.5	Shorter, less complex lesions, PSVR < 2.5	Failed primary efficacy endpoint, PSVR < 2.4	Shorter lesions, fewer CTOs, only 36% of stents deployed nominally, PSVR < 2.0	Shorter, less complex lesions, PSVR < 2.5

<sup>†</sup>For IDE clinical trials with 2 years or more published data as of August, 2020. Results are not directly comparable. Primary patency rates may be defined differently. Information provided is for illustration purposes only, and may differ in head-to-head comparison.

<sup>2</sup>Higher PSVR thresholds may yield higher patency rates. Higashimori A, Kawarada O, Morioka N, et al. Impact of changing PSVR thresholds on the patency rates of SFA recanalisation with self-expanding nitinol stents. *EuroIntervention*. 2013;9(8):964-967.

<sup>1</sup>Matsumura JS, Yamanouchi D, Goldstein JA, et al. The United States Study for Evaluating Endovascular Treatments of Lesions in the Superficial Femoral Artery and Proximal Popliteal by using the Protégé EverFlex Nitinol Stent System II (DURABILITY II). *J Vasc Surg*. July 2013;58(1):73-83.e1.

<sup>2</sup>Rocha-Singh KJ, Bosiers M, Schultz G, et al. A single stent strategy in patients with lifestyle limiting claudication: 3-year results from the Durability II trial. *Catheter Cardiovasc Interv*. July 2015;86(1):164-170.

<sup>3</sup>Laird JR, Katzen BT, Scheinert D, et al. Nitinol stent implantation versus balloon angioplasty for lesions in the superficial femoral artery and proximal popliteal artery: twelve-month results from the RESILIENT randomized trial. *Circ Cardiovasc Interv*. 2010;3(3):267-276.

<sup>4</sup>Laird JR, Katzen BT, Scheinert D, et al. Nitinol stent implantation vs. balloon angioplasty for lesions in the superficial femoral and proximal popliteal arteries of patients with claudication: three-year follow-up from the RESILIENT randomized trial. *J Endovasc Ther*. February 2012;19(1):1-9.

<sup>5</sup>Gray WA, Feiring A, Cioppi M, et al. S.M.A.R.T. self-expanding nitinol stent for the treatment of atherosclerotic lesions in the superficial femoral artery (STROLL): 1-year outcomes. *J Vasc Interv Radiol*. January 2015;26(1):21-28.

<sup>6</sup>Bunte MC, Cohen DJ, Jaff MR, et al. Long-term clinical and quality of life outcomes after stenting of femoropopliteal artery stenosis: 3-year results from the STROLL study. *Catheter Cardiovasc Interv*. 2018;92(1):106-114.

<sup>7</sup>Powell RJ, Jaff MR, Schroë H, Benko A, Diaz-Cartelle J, Müller-Hülsbeck S. Stent placement in the superficial femoral and proximal popliteal arteries with the Innova self-expanding bare metal stent system. *Catheter Cardiovasc Interv*. May 2017;89(6):1069-1077.

<sup>8</sup>Garcia LA, Rosenfield KR, Metzger CD, et al. SUPERB final 3-year outcomes using interwoven nitinol biomimetic supra stent. *Catheter Cardiovasc Interv*. June 1, 2017;89(7):1259-1267.

<sup>9</sup>Ohki T, Angle JF, Yokoi H, et al. One-year outcomes of the U.S. and Japanese regulatory trial of the Misago stent for treatment of superficial femoral artery disease (OSPREY study). *J Vasc Surg*. February 2016;63(2):370-376.e1.

<sup>10</sup>Angle JF, Gasparetto A, Yokoi H, et al. Three-Year Efficacy and Safety of the Misago Peripheral Stent for Superficial Femoral Artery Disease: Final Results from the OSPREY Trial. *J Vasc Interv Radiol*. June 2020;31(6):978-985.

**Indication:** The EverFlex™ Self-expanding Peripheral Stent System is intended to improve luminal diameter in the treatment of symptomatic de-novo or restenotic lesions up to 180 mm in length in the native Superficial Femoral Artery (SFA) and/or proximal popliteal arteries with reference vessel diameters ranging from 4.5–7.5 mm. **Contraindications:** Use of the EverFlex™ Self-expanding Peripheral Stent System is contraindicated in patients with known hypersensitivity to nickel titanium; patients contraindicated for anticoagulant and/or antiplatelet therapy; patients who have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system. **Potential Adverse Events:** Potential adverse events which may be associated with the use of a stent in the SFA and proximal popliteal arteries include, but are not limited to: allergic reaction, amputation, arterial dissection/perforation, bleeding disorders (including GI, lymphatic), infection (local or systemic including bacteremia or septicemia), pseudoaneurysm, restenosis, stent/vessel thrombosis, surgical revascularization.

See the Instructions for Use provided with the product for a complete list of warnings, precautions, adverse events, and device information.

**Caution:** Federal (USA) law restricts these devices to sale by or on the order of a physician.

UC202107072 EN ©2020 Medtronic. All rights reserved. Medtronic and the Medtronic logo are trademarks of Medtronic.™™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 10/2020

**Medtronic**