

# Micra™

## Introducer Sheath

### MI2355A

with hydrophilic coating

## Product specifications

### Introducer Sheath with Hydrophilic Coating

Sheath ID	7.8 mm (23 F)
Working length	55.7 cm (21.9 in)
H x W x D	73 mm x 51 mm x 13 mm

### Dilator

Working length	69.9 cm (27.5 in)
Guidewire compatibility	0.89 mm (0.035 in)

### Additional Features

Stopcock for aspirating and flushing

Radiopaque marker on end of Introducer for location identification



# Medtronic

Medtronic and the Medtronic logo are trademarks of Medtronic.  
™Third party brands are trademarks of their respective owners.  
All other brands are trademarks of a Medtronic company.

**Medtronic**  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA

Toll-free in USA: 800.633.8766  
Worldwide: +1.763.514.4000

**[medtronic.com](https://www.medtronic.com)**

UC201403026b EN ©2020 Medtronic.  
Minneapolis, MN. All Rights Reserved.  
Printed in USA. 02/2020

**Medtronic**