



**Product codes**

**Endurant™ II stent bifurcations**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ETBF         | 23                           | 13                         | C             | 124                       | EE              | 18                           |
| ETBF         | 23                           | 13                         | C             | 145                       | EE              | 18                           |
| ETBF         | 23                           | 13                         | C             | 166                       | EE              | 18                           |
| ETBF         | 23                           | 16                         | C             | 124                       | EE              | 18                           |
| ETBF         | 23                           | 16                         | C             | 145                       | EE              | 18                           |
| ETBF         | 23                           | 16                         | C             | 166                       | EE              | 18                           |
| ETBF         | 25                           | 13                         | C             | 124                       | EE              | 18                           |
| ETBF         | 25                           | 13                         | C             | 145                       | EE              | 18                           |
| ETBF         | 25                           | 13                         | C             | 166                       | EE              | 18                           |
| ETBF         | 25                           | 16                         | C             | 124                       | EE              | 18                           |
| ETBF         | 25                           | 16                         | C             | 145                       | EE              | 18                           |
| ETBF         | 25                           | 16                         | C             | 166                       | EE              | 18                           |
| ETBF         | 28                           | 13                         | C             | 124                       | EE              | 18                           |
| ETBF         | 28                           | 13                         | C             | 145                       | EE              | 18                           |
| ETBF         | 28                           | 13                         | C             | 166                       | EE              | 18                           |
| ETBF         | 28                           | 16                         | C             | 124                       | EE              | 18                           |
| ETBF         | 28                           | 16                         | C             | 145                       | EE              | 18                           |
| ETBF         | 28                           | 16                         | C             | 166                       | EE              | 18                           |
| ETBF         | 28                           | 20                         | C             | 124                       | EE              | 18                           |
| ETBF         | 28                           | 20                         | C             | 145                       | EE              | 18                           |
| ETBF         | 28                           | 20                         | C             | 166                       | EE              | 18                           |
| ETBF         | 32                           | 16                         | C             | 124                       | EE              | 20                           |
| ETBF         | 32                           | 16                         | C             | 145                       | EE              | 20                           |
| ETBF         | 32                           | 16                         | C             | 166                       | EE              | 20                           |
| ETBF         | 32                           | 20                         | C             | 124                       | EE              | 20                           |
| ETBF         | 32                           | 20                         | C             | 145                       | EE              | 20                           |
| ETBF         | 32                           | 20                         | C             | 166                       | EE              | 20                           |
| ETBF         | 36                           | 16                         | C             | 145                       | EE              | 20                           |
| ETBF         | 36                           | 16                         | C             | 166                       | EE              | 20                           |
| ETBF         | 36                           | 20                         | C             | 145                       | EE              | 20                           |
| ETBF         | 36                           | 20                         | C             | 166                       | EE              | 20                           |

**Endurant™ IIIs bifurcations**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ESBF         | 23                           | 14                         | C             | 103                       | EE              | 18                           |
| ESBF         | 25                           | 14                         | C             | 103                       | EE              | 18                           |
| ESBF         | 28                           | 14                         | C             | 103                       | EE              | 18                           |
| ESBF         | 32                           | 14                         | C             | 103                       | EE              | 20                           |
| ESBF         | 36                           | 14                         | C             | 103                       | EE              | 20                           |

**Limbs\***

| Product code |                              |                            |               |                           |                 |                              | Product code   |  |  |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|--|--|--|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) | Total contralateral covered length with EII/EIIs bifurcated† | Total ipsilateral covered length with EIs bifurcated** |  |
| ETLW         | 16                           | 10                         | C             | 82                        | EE              | 14                           | 136  | 155  |  |
| ETLW         | 16                           | 10                         | C             | 93                        | EE              | 14                           | 147  | 166  |  |
| ETLW         | 16                           | 10                         | C             | 124                       | EE              | 14                           | 178  | 177-197  |  |
| ETLW         | 16                           | 10                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 10                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |
| ETLW         | 16                           | 13                         | C             | 82                        | EE              | 14                           | 136  | 155  |  |
| ETLW         | 16                           | 13                         | C             | 93                        | EE              | 14                           | 147  | 166  |  |
| ETLW         | 16                           | 13                         | C             | 124                       | EE              | 14                           | 178  | 177-197  |  |
| ETLW         | 16                           | 13                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 13                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |
| ETLW         | 16                           | 16                         | C             | 82                        | EE              | 14                           | 136  | 135-155  |  |
| ETLW         | 16                           | 16                         | C             | 93                        | EE              | 14                           | 147  | 146-166  |  |
| ETLW         | 16                           | 16                         | C             | 124                       | EE              | 14                           | 178  | 177-197  |  |
| ETLW         | 16                           | 16                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 16                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |
| ETLW         | 16                           | 20                         | C             | 82                        | EE              | 16                           | 136  | 155  |  |
| ETLW         | 16                           | 20                         | C             | 93                        | EE              | 16                           | 147  | 166  |  |
| ETLW         | 16                           | 20                         | C             | 124                       | EE              | 16                           | 178  | 177-197  |  |
| ETLW         | 16                           | 20                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 20                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |
| ETLW         | 16                           | 24                         | C             | 82                        | EE              | 16                           | 136  | 155  |  |
| ETLW         | 16                           | 24                         | C             | 93                        | EE              | 16                           | 147  | 166  |  |
| ETLW         | 16                           | 24                         | C             | 124                       | EE              | 16                           | 178  | 177-197  |  |
| ETLW         | 16                           | 24                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 24                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |
| ETLW         | 16                           | 28                         | C             | 82                        | EE              | 16                           | 136  | 155  |  |
| ETLW         | 16                           | 28                         | C             | 93                        | EE              | 16                           | 147  | 166  |  |
| ETLW         | 16                           | 28                         | C             | 124                       | EE              | 16                           | 178  | 177-197  |  |
| ETLW         | 16                           | 28                         | C             | 156                       | EE              | 16                           | 210  | 209-229  |  |
| ETLW         | 16                           | 28                         | C             | 199                       | EE              | 16                           | 253  | 252-272  |  |

\*The limb mates with the AUI stent graft on the ipsilateral side.

†These calculations assume the minimum 30 mm overlap between the bifurcated stent graft and the contralateral iliac limb per the Endurant II stent graft system *instructions for use*. When using the 124 mm length bifurcated stent graft, subtract 10 mm from total contralateral covered length with bifurcated.

\*\*The 3-5 stent overlap is available only with select limbs. Please refer to the *instructions for use* for more information.

See the device manual for detailed information regarding the *instructions for use*, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.com.

**Aortic extensions**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ETCF         | 23                           | 23                         | C             | 49                        | EE              | 18                           |
| ETCF         | 25                           | 25                         | C             | 49                        | EE              | 18                           |
| ETCF         | 28                           | 28                         | C             | 49                        | EE              | 18                           |
| ETCF         | 32                           | 32                         | C             | 49                        | EE              | 20                           |
| ETCF         | 36                           | 36                         | C             | 49                        | EE              | 20                           |

**Abdominal tubes**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ETTF         | 23                           | 23                         | C             | 70                        | EE              | 18                           |
| ETTF         | 25                           | 25                         | C             | 70                        | EE              | 18                           |
| ETTF         | 28                           | 28                         | C             | 70                        | EE              | 18                           |
| ETTF         | 32                           | 32                         | C             | 70                        | EE              | 20                           |
| ETTF         | 36                           | 36                         | C             | 70                        | EE              | 20                           |

**AUI**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ETUF         | 23                           | 14                         | C             | 102                       | EE              | 18                           |
| ETUF         | 25                           | 14                         | C             | 102                       | EE              | 18                           |
| ETUF         | 28                           | 14                         | C             | 102                       | EE              | 18                           |
| ETUF         | 32                           | 14                         | C             | 102                       | EE              | 20                           |
| ETUF         | 36                           | 14                         | C             | 102                       | EE              | 20                           |

**Iliac extensions**

| Product code |                              |                            |               |                           |                 |                              |
|--------------|------------------------------|----------------------------|---------------|---------------------------|-----------------|------------------------------|
|              | Proximal graft diameter (mm) | Distal graft diameter (mm) | Distal design | Total covered length (mm) | Delivery system | Catheter outer diameter (Fr) |
| ETEW         | 10                           | 10                         | C             | 82                        | EE              | 14                           |
| ETEW         | 13                           | 13                         | C             | 82                        | EE              | 14                           |
| ETEW         | 20                           | 20                         | C             | 82                        | EE              | 16                           |
| ETEW         | 24                           | 24                         | C             | 82                        | EE              | 16                           |
| ETEW         | 28                           | 28                         | C             | 82                        | EE              | 18                           |

See the device manual for detailed information regarding the *instructions for use*, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at [medtronic.com](http://medtronic.com).

For applicable products, consult *instructions for use* on [www.medtronic.com/manuals](http://www.medtronic.com/manuals). Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Important Reminder: This information is intended only for users in markets where Medtronic products and therapies are approved or available for use as indicated within the respective product manuals. Content on specific Medtronic products and therapies is not intended for users in markets that do not have authorization for use.

# Medtronic

[Medtronic.com/aortic](http://Medtronic.com/aortic)

UC202301796 ML ©2022 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. Not for distribution in the USA or France. 08/2022