

PRODUCT SPECIFICATIONS

Evolut™ PRO+ TAVR

BIOPROSTHESIS



Evolut PRO+
23 mm Valve

Model Number:
EVPROPLUS-23US

Size: 23 mm

Aortic annulus
diameter:
17*/18-20 mm



Evolut PRO+
26 mm Valve

Model Number:
EVPROPLUS-26US

Size: 26 mm

Aortic annulus
diameter:
20-23 mm



Evolut PRO+
29 mm Valve

Model Number:
EVPROPLUS-29US

Size: 29 mm

Aortic annulus
diameter:
23-26 mm



Evolut PRO+
34 mm Valve

Model Number:
EVPROPLUS-34US

Size: 34 mm

Aortic annulus
diameter:
26-30 mm

LOADING SYSTEM



Evolut PRO+
Loading System

Model Number:
L-EVPROP2329US

23, 26, and 29 mm TAV



Evolut PRO+
Loading System

Model Number:
L-EVPROP34US

34 mm TAV

*Measurement for TAV-in-SAV only.
All items above are single-use. Refer to Instructions for Use for a complete list of warnings, precautions, and contraindications.

DELIVERY CATHETER SYSTEM

Evolut PRO+


Delivery Catheter System with InLine Sheath

14 Fr Equivalent (6.0 mm Capsule Outer Diameter)

Model Number:

D-EVPROP2329US

23, 26, and 29 mm TAV



Evolut PRO+


Delivery Catheter System with InLine Sheath

18 Fr Equivalent (7.33 mm Capsule Outer Diameter)

Model Number:

D-EVPROP34US

34 mm TAV



Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel: (763) 514-4000
Fax: (763) 514-4879
Toll-free: (800) 328-2518

LifeLine
CardioVascular Technical Support
Tel: (877) 526-7890
Tel: (763) 526-7890
Fax: (763) 526-7888
rs.cstechsupport@medtronic.com

©2019 Medtronic.
All rights reserved.
Medtronic, Medtronic logo,
and Further, Together are
trademarks of Medtronic.
All other brands are trademarks
of a Medtronic company.

[medtronic.com/TAVR](https://www.medtronic.com/TAVR)

UC202001772 EN
09/2019